


# Pediatric Wheelchair Parts Guide

---



## Description of Components

**Push Handle:** Used by a caregiver to push or maneuver the wheelchair.

**Headrest:** Provides support and stability for the user's head and neck.

**Backrest:** Supports the user's back and helps maintain posture.

**Armrest:** Provides support for the user's arms and shoulders, aiding in comfort and weight distribution.

**Seat Cushion:** Provides a comfortable base for the user to sit on and helps prevent pressure sores.

**Wheel Lock (Brake):** Used to lock the rear wheels in place to prevent the chair from rolling when transferring or stationary.

**Footrest Hanger:** Supports the footplate and can be adjusted for height and angle.

**Footplate:** The platform where the user places their feet.

**Front Caster Wheel:** The smaller wheels at the front that allow for easy steering and turning.

**Rear Wheel (Tire):** The larger main wheels that provide movement and stability.

**Anti-Tip Tube:** A small safety feature at the back that prevents the wheelchair from tipping backward.

**Disclaimer:** This list is not all-inclusive, and not all wheelchairs are the same. Features, configurations, and components may vary significantly based on the specific model, custom requirements, and intended use. Always consult with a qualified professional, such as a physical or occupational therapist, for individualized equipment recommendations and setup.